

2× qPCR Mix for microRNA (ROX1 plus)

Catalog No.: EZB-miQP1

Description

The EZBioscience® 2× qPCR Mix for microRNA (ROX1 plus) uses a specially modified *Taq* DNA polymerase protected by a hot-start activation technique, and optimized qPCR buffer system to perform SYBR Green I based Real-time quantitative PCR (qPCR) for microRNA expression detection. The mix is prepared at 2× reaction concentration, and is specially optimized for microRNA, which strikingly improves the microRNA amplification and reduces the non-specific amplification. Meanwhile, the ROX1 passive reference dye (high ROX) is premixed in the mixture to make it an excellent ready to use qPCR reagent for microRNA expression detection. Besides, a universal 3' qPCR primer is also supplied in this Kit. It is recommended that the kit be used together with the EZBioscience® EZ-press Cell to cDNA Kit PLUS for microRNA, or microRNA Reverse Transcription Kit (Cat. No.: EZB-miRT1 or EZB-miRT2).

Components

| Components | EZB-miQP1-S (200 Rxns) | EZB-miQP1-L (500 Rxns) |
|----------------------------------|--------------------------|--------------------------|
| 2× qPCR Mix for microRNA* | 2 ml (1 ml × 2 tubes) | 5 ml (1 ml × 5 tubes) |
| Universal 3' qPCR Primer (10 μM) | 100 μl | 250 μl |
| U6 Primer (10 μM) | 50 μl | 125 μl |

*: Contain Hot-start DNA Polymerase, dNTPs, Mg²⁺, SYBR Green I dye, and premixed with and ROX1, to rectify the error of fluorescence signals between different wells.

Storage

All components should be stored at -20°C and avoid of light. The quality of this product is guaranteed for 18 months under recommended storage conditions.

Notice

Appropriate Real-time PCR instruments for this Kit are as follows:

| |
|--|
| ABI 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast; StepOne™, StepOne Plus™. |
| Bio-Rad CFX96™, CFX384™, iCycler iQ™, iQ™5, MyiQ™, MiniOpticon™, Opticon®, Opticon 2, Chromo4™; Cepheid SmartCycler®; Roche LightCycler™ 480; Eppendorf Mastercycler® ep realplex, realplex 2 s; Illumina Eco qPCR; Qiagen/Corbett Rotor-Gene® Q, Rotor-Gene® 3000, Rotor-Gene® 6000; Thermo Scientific PikoReal Cycler. |

Protocol (Suggestions for qPCR Amplification)

1. The suggested template amount is 4 μl cDNA template (the cDNA derived from the EZBioscience® microRNA Reverse Transcription Kit, Cat. No.: EZB-miRT1 or EZB-miRT2), then diluted for 5 times, which is 100 μl. The following loading table is designed for reaction volume of both 10 μl and 20 μl:

| Components | Amount/ reaction (μl) | Amount/ reaction (μl) | Final concentration |
|--------------------------------------|-----------------------|-----------------------|---------------------|
| Total Reaction Volume | 10 | 20 | - |
| 2× qPCR Mix for microRNA | 5 | 10 | 1× |
| Specific forward qPCR Primer (10 μM) | 0.2 | 0.4 | 100 nM |
| Universal 3' qPCR Primer (10 μM) | 0.2 | 0.4 | 100 nM |
| Template cDNA | 2 | 4 | 1 ~ 100 ng |
| ddH ₂ O | 2.6 | 5.2 | - |

2. Standard 2-step amplification program:

| Step | 1 | 2 | |
|--------|-------------------------------------|-------------------|----------------|
| | Hot-Start DNA Polymerase Activation | PCR | |
| | Hold | Cycle (40 cycles) | |
| | | Denature | Anneal/ Extend |
| Temp | 95°C | 95°C | 60°C |
| Time | 5 mins | 10 secs | 30 secs |
| Volume | 50 μl/ 20 μl | | |

Melt Curve could be set following the recommendation of the equipment. A representative program for Melt Curve could be set up as follows:

| Step | 1 | 2 | 3 |
|-----------------------|--------|-------|---------------------------------|
| Heating/ Cooling rate | 100% | 100% | 1% |
| Temp | 95°C | 60°C | 95°C |
| Time | 15 sec | 1 min | 30 sec |
| Data Collection | - | - | During temperature rising stage |

Important: The 95°C for 5 min step is required to activate the Hot-Start DNA Polymerase.

Data Analysis

Two types of quantification methods are commonly applied to quantify target microRNA expression when using this SYBR Green qPCR Mix.

1. Relative Quantitation: target microRNA expression is measured against an internal standard, such as U6.

microRNA expression can be measured by the quantitation of cDNA converted from a microRNA corresponding to this microRNA relative to a calibrator sample serving as a physiological reference. In a typical experiment, microRNA expression levels are studied as a function of either a treatment of cells in culture, of patients, or of tissue type, or of exosomes. The calibrator sample in each case is the cDNA from either the untreated cells or patients, or a specific tissue type, or of exosomes. All quantitations are also normalized to an endogenous control such as U6 to account for variability in the initial concentration and quality of the total RNA and in the conversion efficiency of the reverse transcription reaction.

Set Ct^{A1} as the Ct value of the target microRNA of sample 1, and Ct^{B1} as the Ct value of the internal control gene of sample 1; set Ct^{A2} as the Ct value of the target microRNA of sample 2, and Ct^{B2} as the Ct value of the internal control gene of sample 2. The expression difference (in folder) of the target microRNA in sample 1 and in sample 2 can be calculated this way (2^{-ΔΔCt} approach): ΔΔCt = (Ct^{A2} - Ct^{B2}) - (Ct^{A1} - Ct^{B1})

The expression level of the target microRNA in sample 2 is 2^{-ΔΔCt} times that of sample 1.

Note: This calculation method is based on the assumption that the amplification efficiency is 100% (the amount of products after each cycle is doubled). If the amplification efficiency is not 100%, the calculation formula needs to be

amended according to actual reaction efficiency. For example: if the amplification efficiency of the target microRNA and internal control genes is 1.90, then the formula should be amended into $1.90^{\Delta\Delta C_t}$.

2. Absolute quantitation: compares the Ct of an unknown sample against a standard curve with known copy numbers.

Absolute quantitation is applicable only if isolation procedure and sample contents do not affect PCR amplification.

Determination of baseline and threshold: please refer to the real-time PCR system software used to calculate baseline and threshold values for a detector or manually set up according to amplification curves.

Passive reference ROX: is a dye molecule included in the SYBR Green PCR Master Mix that does not participate in the PCR amplification. On applied Biosystems real-time PCR system, the passive reference provides an internal reference to which the SYBR Green/ dsDNA complex signal can be normalized during data analysis.

Representative Results




Figure 1. miR-21, miR-1 and U6 expression level in cDNA from A549 cells were detected. The amplification plots and dissociation curve above indicate that, all the above microRNAs could be amplified efficiently and specifically. (Notes: Quantitation is fine when Ct value is within 13 ~ 32; if the Ct value is larger than 35, the amplification is invalid; If the Ct value is too small (<13), please dilute the template and repeat the quantitation; if the Ct value is within 32 ~ 35, 3 replicates are required to validate the results. A typical valid dissociation curve should be unimodal due to highly specific amplification and no primer dimer formation (Figure 1). If the dissociation curve is apparently multimodal, it often suggests primer dimer formation, non-specific amplification, or other contamination. Alternatively, confirmation of non-specific amplification can be analyzed by agarose gel electrophoresis).

Suggestions for Primer Design

The reverse primer for microRNA quantification (a universal qPCR primer) is already provided in this Kit, which is corresponding with the EZBioscience® microRNA Reverse Transcription Kit (Cat. No.: EZB-miRT1 or EZB-miRT2). The

following general principals may be considered when designing primers for target microRNAs:

1) First, get the target microRNA sequence from microRNA data base (such as miRBase: <http://www.mirbase.org/>).

2) Copy the target microRNA sequence, change the U in the sequence to T, then delete the last 6 nucleotides at the 3' end.

3) Add 3 ~ 6 nucleotides to the 5' end of the primer, to justify the Tm value of the primer to 55 ~ 60°C (the added sequence mainly contains G&C, such as CGGGC, GCGGGC, or A/TGCCCG).

4) Synthesize the primer and exam the quality of the primer by qPCR. If the Ct value derived is between 15 ~ 32, the amplification curve displays representative S form curve, and the melt curve has a typical single peak, the primer could be approved. If there are apparently double peaks in the melt curve, optimize the primer by add or delete the specific nucleotides to the 5' or 3' end of the primer.

Trouble Shooting

| Problem | Potential Cause(s) | Suggestion(s) |
|---|---|--|
| Apparent amplification can be observed in negative control. | The reagents or water used is contaminated. | Change new reagents or water and retry. The reaction should be set up in a super clean bench to minimize contamination from the air. |
| | Appearance of primer dimer. | It is normal that amplification occurs in negative control after 35 cycles. Analysis can be performed combining with the melt curve. |
| Ct value appears too late (large). | Low amplification efficiency. | Optimize the reaction. Try three-step program or re-design the primers. |
| | The concentration of templates is too low. | Increase the amount of the template, meanwhile avoid destroying the amplification. |
| | Degradation of templates. | Prepare new templates and retry. |
| Abnormal shape of amplification plot. | There are PCR inhibitors in the reaction. | Inhibitors are usually brought in when adding templates. Increase the dilution folds or prepare new templates and retry. |
| | Rough amplification plot. | It is caused by system rectification due to weak signal. Elevate the template concentration and repeat the reaction. |
| | Broken or downward amplification plot. | The concentration of templates is too high. The end value of the baseline is bigger than Ct value. Decrease the end of the baseline (Ct value - 4), and re-analyze the data. |
| No amplification plot. | Amplification plot goes downward suddenly. | There are bubbles left in the reaction tube, which break up when the temperature rises, thus the instrument detects the sudden decrease of the fluorescence value. Spin briefly and check closely if there are bubbles left before reaction. |
| | Cycling number is insufficient. | Generally, the cycling number is set to be 40. But notice that too many cycles will result in excessive background, thus reducing the data reliability. |
| | Check if there is signal collection procedure during cycling. | In two-step program, signal collection is usually positioned at annealing and extension stage. |
| | Check if the primers are degraded. | Test the integrity of primers after long-term storage through PAGE electrophoresis to confirm if the primers are degraded already. |
| | The concentration of templates is too low. | Reduce the dilution fold and retry. For target gene with unknown expression level, begin without template dilution. |
| Degradation of templates. | Prepare new templates and retry. | |